


---

# **Real-time PCR (qPCR)**

Dr. Luu Phuc Loi, PhD.  
Presenter: Hoang Kim


# The Evolution: From End-Point to Real-Time Detection

## The Past: Conventional PCR



- **Method:** End-point detection (Agarose Gel)
- **Nature:** Qualitative / Semi-quantitative
- **Drawbacks:** Labor-intensive, Open-tube system (Contamination risk)
- **Insight:** Analysis of the 'remains' after the reaction.

## The Present: Real-Time PCR (qPCR)



- **Method:** Kinetic detection (Fluorescence)
- **Nature:** Quantitative, High-throughput
- **Advantages:** Closed-tube system, Wide dynamic range
- **Insight:** Observation of the reaction 'live' as it happens.

# The Problem with End-Point Analysis: The Plateau Effect

## The Phenomenon


As reagents (dNTPs, primers) deplete and enzyme activity slows, the reaction stops generating template exponentially.

## The Flaw


End-point band intensity correlates with *yield*, not the *initial input*. Looking at Cycle 40 (the plateau), Sample 1 and Sample 3 appear identical despite different starting quantities.

## The Analogy


You cannot judge the speed of runners by taking a photo at the finish line after everyone has stopped running.



# Decoding the Kinetics: The Four Phases of PCR



# The Exponential Phase: The Window of Quantification



## Why do we measure here?

### 1. Non-Limiting Reagents

dNTPs and primers are in excess. The reaction is not starving.

### 2. Peak Efficiency


The polymerase is working at maximum capacity.

$$\text{Equation: } \text{Amount} = \text{Start} \times 2^{\text{cycles}}$$


### 3. Precision

The relationship between Fluorescence ( $R_n$ ) and Template Amount is linear on a log scale ( $R^2 > 0.99$ ). We ignore the plateau and the baseline.

# Setting the Parameters: Baseline and Threshold



# The Fundamental Metric: The Ct Value



## Definition

The fractional cycle number at which the fluorescence signal crosses the threshold line.


## Synonyms

- $C_q$  (Quantification Cycle) – MIQE preferred
- $C_p$  (Crossing Point) – Roche LightCycler


## Key Concept

$C_t$  is a time-based measurement (when did it happen?), not an intensity-based measurement.

# The Inverse Relationship



High Template Load



Low Ct Value


Low Template Load



High Ct Value

**Rule: Ct is inversely proportional to the log of the initial template copy number.**  
 $\Delta 3.3 \text{ cycles} \approx 10\text{-fold difference.}$

# Dynamic Range and Sensitivity



## Sensitivity:

Detection down to 1 copy.


## Precision:

Reliable differentiation between 10,000 and 20,000 copies.



## Linearity:

Linear results across 8 orders of magnitude.

# Factors Influencing the $C_t$ Value



# Normalization: The Role of Passive Reference Dyes (ROX)



## The Problem:

Pipetting errors, air bubbles, and optical variations cause signal fluctuations.

## The Solution:

Passive Reference Dye (ROX). A dye that does not participate in PCR provides a constant internal control.

## The Formula:

$$R_n = \frac{\text{Emission of Reporter Dye (FAM)}}{\text{Emission of Passive Reference (ROX)}}$$

# Module 1 Summary: The Theoretical Pillars



## Quantification Window

Analysis must occur in the **Exponential Phase**, where reagents are non-limiting and efficiency is constant.




## The Metric

The  $C_t$  (or  $C_q$ ) **value** is the gold standard for quantification, representing the cycle where signal exceeds noise.



## The Relationship

$C_t$  is **inversely proportional** to the starting template quantity. ( $\$$ High Load = Low  $C_t$ ).



## Sensitivity

qPCR offers a dynamic range of 7–8 logs, vastly superior to end-point gel analysis.



## Normalization


Reference dyes (ROX) and baseline correction are essential for removing non-biological noise.

# Two Strategies to Visualize the Invisible

## Classification of Detection Chemistries

### 1. Non-Specific Detection


(The “Floodlight” Strategy)



- **Mechanism:** DNA-binding dyes intercalate into any double-stranded DNA (dsDNA).
- **Principle:** ‘If it is double-stranded, it glows.’
- **Key Example:** SYBR® Green I

### 2. Specific Detection


(The “Sniper” Strategy)



- **Mechanism:** Sequence-specific oligonucleotide probes labeled with fluorophores.
- **Principle:** ‘Glow only when the specific target code is found.’
- **Key Examples:** Hydrolysis Probes (TaqMan®), Molecular Beacons

# The Mechanism of SYBR® Green I

From Darkness to Light: Intercalation



**The Correlation:** Fluorescence intensity is directly proportional to the amount of dsDNA amplicon accumulated.

# Why Choose SYBR Green?

## The Advantages of DNA-Binding Dyes



### Cost-Effective

Significantly less expensive than labeled probes. Ideal for large-scale screening or limited budgets.



### Simplicity

Requires only standard PCR primers. No need to design complex internal probes or labeled oligos.



### Universality

"One dye fits all." The same master mix works for any gene target, from GAPDH to viral pathogens.

**\*\*Ideal Use Case:\*\* Initial expression validation, primer optimization, and high-throughput screening.**

# The ‘Blindness’ of SYBR Green

## Limitations: Non-Specificity and Artifacts

### The Problem:

SYBR Green binds indiscriminately to *any* double-stranded DNA. It cannot distinguish between the target gene and amplification artifacts.


### Common Artifacts:

- **Primer-Dimers:** Primers binding to each other.
- **Non-Specific Products:** Mis-priming elsewhere in the genome.

### The Consequence:

False positives and overestimation of target quantity (Ct values appear too early).

### What the Machine Sees vs. Reality



**Crucial Requirement:** Users must verify reaction specificity post-amplification using Melt Curve Analysis.


# The Solution: Melt Curve Analysis

Using Temperature to Determine Identity

**The Protocol:** Performed immediately after PCR cycles.

1. **Ramp:** Temperature is slowly increased ( $60^{\circ}\text{C} \rightarrow 95^{\circ}\text{C}$ ).
2. **Denaturation:** dsDNA destabilizes and separates into single strands.
3. **Release:** Dye is released back into solution.
4. **Signal Drop:** Fluorescence decreases rapidly.


**Definition:**  $T_m$  (Melting Temperature) is the point where 50% of the amplicon is denatured.



# Interpreting the Dissociation Curve


Identifying Success vs. Contamination

**Scenario A: Clean Reaction (Pass)**



**Specific Amplification.** Single peak indicates one pure product.


**Scenario B: Contamination (Fail)**



**Non-Specific Binding / Primer-Dimers.** Multiple peaks or shoulders indicate artifacts.

# Hydrolysis Probes (TaqMan®) Structure


## The Components of Specificity



**Design:** The probe binds specifically to the target DNA between the forward and reverse primers.

# The Mechanism: FRET

## Fluorescence Resonance Energy Transfer




**The Setup:** During annealing (55-65°C), the probe binds specifically to the target. Because the probe is intact, the Quencher suppresses the Reporter. The signal is silent.

# Signal Generation: The Role of Taq Polymerase


$5' \rightarrow 3'$  Exonuclease Activity

## Step 1: Encounter




Polymerase Extends

## Step 2: Hydrolysis



$5' \rightarrow 3'$  Exonuclease Activity


## Step 3: Signal Release



FRET Broken = **Fluorescence**

# Multiplexing Capability

Multiple Targets, One Tube



**The Power of Color:** Distinct probes with different reporter dyes allow simultaneous detection.

Benefits:

- **Internal Controls:** Normalize data within the same well.
- **Efficiency:** Conserves sample and reagents.
- **Throughput:** Detect multiple pathogens at once.


# Comparison: SYBR® Green vs. TaqMan®

## Choosing the Right Chemistry

Feature	SYBR® Green (Dye)	TaqMan® (Probe)
Specificity	Low (Primers only)	High (Primers + Probe)
Cost	Low (Entry-level friendly)	High (Requires custom probe)
Assay Design	Simple (Primers only)	Complex (Probe design needed)
Multiplexing	No	Yes (Multicolor)
Post-PCR Analysis	Required (Melt Curve)	Not Required


# Defining the Goal: Absolute vs. Relative Quantification

## Absolute Quantification




Question: "How many copies are there?"

- Determines exact copy number or unit mass.
- Requires interpolation against a Standard Curve.
- Example: Viral Load (Copies/mL).



## Relative Quantification




Question: "How much did it change?"

- Calculates fold-change ratios.
- Requires normalization to Reference Gene & Calibrator.
- Example: Gene Expression (5-fold increase).

# Absolute Quantification: Calibrating the Ruler

Measuring the Unknown against the Known




**Standard Templates:**  
Plasmid with cloned GOI  
Synthetic Oligonucleotides  
Genomic DNA  
In Vitro Transcripts

**Critical Rule:** Quantification must be by Interpolation (within the curve), never Extrapolation.

# Constructing the Standard Curve

## The Serial Dilution Workflow




### Design Requirements:


- **Dynamic Range:** Must cover the expected range of unknowns.
- **Points:** Minimum 4 concentration points.
- **Replicates:** Duplicates or Triplicates to reduce pipetting error.

# Quality Control: Evaluating Linearity ( $R^2$ )

Is the ruler straight?



Excellent Linearity




Poor Fit / Non-Linear

**$R^2$  (Pearson Correlation Coefficient):**  
Measures fit quality.  
**Acceptance Criteria:**  $> 0.980$   
**⚠ Warning:** Bending at high Ct  
indicates Limit of Detection (LOD).

# PCR Efficiency (E): The Doubling Principle

Theoretical Perfection vs. Biochemical Reality



## The Ideal Slope: -3.32

Calculated from:


$$10\text{-fold increase} = 3.32 \text{ cycles } (2^{3.32} \approx 10)$$

**Acceptable Range:**

90% – 105%

**Slope Range:**

-3.1 to -3.6



10 fold change =  $2^{3.3}$  i.e. a 3.3 cycle difference.

PCR double each cycle =  $2^X$  where x is number cycles

# Calculating Efficiency from Slope

Deriving the percentage from the regression line

$$\text{Efficiency } (E) = 10^{-1/\text{slope}} - 1$$

Derived from Standard Curve

## Worked Example


Given Slope =  $-3.32$

1. Step 1: Divide  $-1$  by  $-3.32 = 0.301$
2. Step 2: Calculate  $10^{0.301} = 2.0$
3. Step 3: Subtract  $1 = 1.0$

**Result: 100% Efficiency**


# Diagnosing Efficiency Issues

What the numbers reveal about reaction health



# Relative Quantification: Measuring Gene Expression

Focus: Fold-Change Ratios




## Objective:

- Determine up-regulation or down-regulation.
- Requires a Calibrator (Untreated/Time Zero).
- CRITICAL: Must be normalized to a Reference Gene to account for loading errors.

# The Anchor: Selecting Reference Genes

Normalizing for Error



## The Role:

- Corrects for variations in RNA input.
- Corrects for Reverse Transcription efficiency.

## The Criteria:


- Constitutive expression.
- **STABLE** across all treatments (Drug, Temp, Time).

## Common Anchors:

- GAPDH
- Beta-Actin
- 18S rRNA

**Validation is key.**

If the anchor moves (expression changes), the data is invalid.



# Step 1: The Delta Ct Calculation

Normalizing Within the Sample

$$\Delta\text{Ct} = \text{Ct}(\text{Target Gene}) - \text{Ct}(\text{Reference Gene})$$

- This calculation is performed for **every** sample individually.
- Goal: Eliminate loading bias. If you loaded 2x RNA, both Ct values shift, but the difference ( $\Delta\text{Ct}$ ) remains constant.

Sample	Target Ct	Ref Ct	$\Delta\text{Ct}$
Treated	24.5	18.0	6.5
Control	26.5	18.0	8.5

## Step 2: The Delta-Delta Ct Calculation

Comparing Sample vs. Calibrator

$$\Delta\Delta Ct = \Delta Ct(\text{Sample}) - \Delta Ct(\text{Calibrator})$$



Usually the Untreated Control

- This subtracts the baseline noise of the control from the sample signal.
- If Result is 0: No difference.
- If Result is Negative: Target is Up-regulated (Lower Ct = More Template).

Example Continuation

$$\Delta\Delta Ct = 6.5 \text{ (Treated)} - 8.5 \text{ (Control)} = -2.0$$

## Step 3: The Fold Change Formula (Livak Method)

The  $2^{\Delta\Delta Ct}$  Method

$$\text{Fold Change} = 2^{-\Delta\Delta Ct}$$

Assumption of 100%  
Efficiency (Doubling)

The difference in cycles  
calculated in Step 2.

Example Calculation

$$\begin{aligned}\text{Fold Change} &= 2^{-(-2.0)} \\ &= 2^2 \\ &= \textbf{4-Fold Increase}\end{aligned}$$

### Founders Grotesk

- Efficiencies of Target and Reference must be approx. EQUAL and near 100%.

### Crucial Assumption

- Efficiencies of Target and Reference must be approx. EQUAL and near 100%.
- ⚠ If efficiencies differ, **this method fails**.

# The Pfaffl Method

When Efficiencies Differ

What if Target E = 90% and Reference E = 100%?

The Livak method will be inaccurate.

$$\text{Ratio} = \frac{(E_{\text{target}})^{\Delta Ct(\text{control-sample})}}{(E_{\text{ref}})^{\Delta Ct(\text{control-sample})}}$$


Uses specific Efficiency values calculated from standard curves instead of assuming "2".

## Takeaway


Use Pfaffl for complex biological samples or when perfect PCR conditions are impossible.

# Data Analysis Workflow Summary

From Validation to Calculation



# Primer Design 101: The Foundation of Specificity



Ideally:  $\Delta G$  of Dimer Formation > -4 kcal/mol (Avoid Self-Binding).


# Amplicon Design: Optimization for Efficiency

## OPTIMAL: The Sprinter



- **Target Length:** 70–150 bp (Probe) / 100–300 bp (SYBR)
- **Benefit:** Polymerase completes extension before cycle ends.
- **Condition:** Spans Exon-Exon Junction (cDNA only).


## SUB-OPTIMAL: The Marathon Runner



- **Target Length:** > 300 bp
- **Risk:** Secondary structures (hairpins) block binding.
- **Result:** Reduced Efficiency.

**PRO TIP:** Use 'mfold' software to predict and avoid secondary structures.

# Troubleshooting: The ‘Flatline’ (No Amplification)



## REAGENTS (Chemistry)



Missing Taq, dNTPs, or Primers?

**Action:** Check Master Mix.  
Run Positive Control.

## TEMPLATE (Sample)



Degraded or Low Concentration?

**Action:** Check integrity on Bioanalyzer.


## INHIBITION (Environment)



Ethanol, Salt, or Hemoglobin carryover?


**Action:** Dilute Template (1:10). If Ct improves, inhibition is present.

# Troubleshooting: Delayed Amplification (Late Ct)



# Clinical Application: Absolute Quantification (Viral Load)

## The Standard Curve



## The Patient Report


### Viral Load Assay: Hepatitis B (HBV)

Sample ID: #89204-X

**$1.5 \times 10^5$  copies/mL**



Status: Detected - Monitor Therapy Response

# Diagnostic Gold Standard: COVID-19 Detection



# Genotyping: The Mechanism of SNP Detection

TaqMan Probe Competition (Allelic Discrimination)




Match = Cleavage = Signal

Mismatch = Displacement = No Signal


MGB Probes allow for shorter sequences (~13bp), increasing mismatch sensitivity.

# Data Analysis: Reading the Allelic Discrimination Plot



Endpoint Analysis: Position on the plot determines genotype identity.

# Agricultural Application: GMO Quantification



$$\% \text{ GMO} = \frac{\text{Copy Number of Target Gene}}{\text{Copy Number of Reference Gene}} \times 100$$

(e.g., Lectin - Total Organism DNA)



# Environmental Monitoring: Water & Soil Safety

## Water Quality



**Target:** Fecal Coliforms (*Enterococcus*)

↳ **Benefit:** Speed: Result in hours (vs. 24+ days for culture).

**Context:** Beach & Recreational Safety.

## Bioremediation



**Target:** *alkB* gene (Alkane Hydroxylase)

↳ **Benefit:** Tracking oil-eating bacteria populations.

**Context:** Confirming cleanup efficacy.

---

**THANKS FOR  
LISTENING**

---